

CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175

CERTIFICATE & REPORT NO: INV/EN14175/562-B
ISSUE DATE: 24th April 2013

Fume Cupboard Type:

chemtop1200

CF-120-PP ductless bench type

External Dimensions:

Height = 2060mm
(1260mm + 800mm)

Width = 1200mm

Depth = 700mm (max)

Internal Dimensions:

Height = 700mm (max)

Width = 1125mm

Depth = 620mm (max)

Baffle Gap Dimensions:

Top: 32x75mmx8mm slots

Bottom: 116x75mmx8mm slots

Fume Cupboard Manufacturer:

TopAir USA

3182 Monterey Drive

Merrick

NY 11566

USA



Fume Cupboard Flow:

Volume flow rate: 0.310m³/s

Test opening width: 1125mm

Test opening height: 500mm

Face velocity: 0.520m/s

Pressure drop: N/A

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

Tested and Certified by: Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.