

## CERTIFICATE OF TYPE TESTING IN ACCORDANCE WITH EN 14175

CERTIFICATE & REPORT NO: INV/EN14175/556-A  
ISSUE DATE: 24th April 2013

**Fume Cupboard Type:**

fumetop1500

**FH-150 bench type**

*External Dimensions:*

Height = 2300mm  
(1500mm + 800mm)

Width = 1500mm

Depth = 800mm

*Internal Dimensions:*

Height = 1050mm

Width = 1250mm

Depth = 670mm (wall-sash)  
620mm (baffle-sash)

*Baffle Gap Dimensions:*

Top gap: 40mm

Bottom gap: 140mm - 240mm

**Fume Cupboard Manufacturer:**

TopAir USA

3182 Monterey Drive

Merrick

NY 11566

USA



**Fume Cupboard Flow:**

**Volume flow rate:** 0.345m<sup>3</sup>/s

**Test opening width:** 1250mm

**Test opening height:** 500mm

**Face velocity:** 0.511m/s

**Pressure drop:** 66Pa @ 250mmD

This is to certify that the fume cupboard described above has been type tested in accordance with Part 3, in compliance with the requirements of Part 2 and with reference to Part 1 of EN 14175 and resulted in performance characteristics given in the corresponding test report.

**Tested and Certified by:** Dr A F Bicen

Invent UK Ltd, 85 Southdown Road, Harpenden, Herts AL5 1PR

Tel: 01582 461666 : Fax: 01582 460077 : E-mail: info@invent-uk.com

A handwritten signature in blue ink that reads 'A.F. Bicen'.